

entry	
editor	
process	
MRN	

MEASUREMENT SHEET

partial hand prosthesis V1

order number	
company	
street	
technician	

date	
e-mail	
zip city	
telephone	

patient information


Name	
date of birth	

affected side	left	right
gender	male	female

notice

configuration

- | | | | | |
|---|--------------------------|-------------|----------|-------|
| trail prosthesis silicone | partial hand with finger | fingernails | silicone | acryl |
| partial hand prosthesis basic | complete hand | | | |
| partial hand prosthesis classic (color ring will be sent) | | | | |
| partial hand prosthesis individual (patient have to come to dresden for color matching) | | | | |



- scribble stump in sketch
- measure the circumference of the stump and enter it into the sketch at the measuring points
- measure the circumference of the adjacent fingers and corresponding fingers from the obtained hand and scribble it in sketch

